



Axlab BLADES



Item
1993-01

Plasma
Low profile

50 pcs

Axlab Blades **Universal35**
Low-profile blade.
For most sectioning tasks.
Box with 10 dispensers,
500 blades total.

Item
1993-02

Plasma
Low profile

50 pcs

Axlab Blades **Solid35**
Low-profile blade.
For dense or hard tissue.
Box with 10 dispensers,
500 blades total.

Item
1993-03

Plasma
High profile

50 pcs

Axlab Blades **Versatile35**
High-profile blade. For frozen
sectioning, mega cassettes and
trimming. Box with 10 dispensers,
500 blades total.

High-quality consumables,
crafting precision for efficient
and reliable laboratory workflows.

AXLAB.DK

